

MGI MELTBLOWN HEALTH MASK INSERT

Mask Insert with 5 Layers

Helps filter particles,
droplets, and air pollution.



5 LAYERS

- Non-woven Layer
- Melt blown Layer
- Activated Carbon Layer
- Melt blown Layer
- Non-woven Layer



5 Layers



Replace Every
72 Hours



Practical &
Hygienic



For optimal use, wash your hands before inserting.

MGI MELTBLOWN HEALTH MASK INSERT

5-Layer Insert for Additional Protection




An additional insert placed inside the mask to help filter droplets, dust, and other harmful airborne particles.



Recommended to be replaced every 72 hours and wash your hands before installing the insert.



5-LAYER CONSTRUCTION



BENEFITS



Helps filter droplets and airborne particles



Helps block dust, bacteria, and harmful substances



Comfortable to use inside the mask



Extra protection for daily activities



How to Use:

Insert the insert into the mask before use.



Replacement:

Replace the insert every 72 hours for optimal results.



Note:

Wash your hands thoroughly before installing the insert.



5 LAYERS



Additional Protection



Practical & Hygienic